

# Practical and Technical Advantages of DXA Visceral Fat Assessment Compared with Computed Tomography

Thomas L. Kelly, Senior Principal Scientist  
Hologic, Inc.

## Abstract

Over the last decade obesity has exploded to epidemic proportions in America. Obesity-related health risks include coronary heart disease, type 2 diabetes, stroke, dyslipidemia, and certain types of cancers.

Traditional anthropomorphic measures of obesity like BMI and waist circumference quantify excess weight or size, not excess fat, and misclassify patients by failing to assess pathogenic visceral fat depots. DXA visceral fat measurements may be superior to traditional measures for evaluating obesity-related disease risk. The available literature supports visceral fat thresholds of 100 cm<sup>2</sup> for increased risk and 160 cm<sup>2</sup> for high risk and suggests this classification scheme may identify patients most likely to benefit from preventive interventions.

## Introduction

Obesity is widely recognized as America's most serious public health problem following a Center for Disease Control (CDC) report that 78 million adults (37.5%) and 12.5 million children (17%) are obese.<sup>1</sup> For the first time in recorded history, obesity threatens to reduce life expectancy from both chronic and acute diseases including cardiovascular disease, stroke, dyslipidemia, cancer, hypertension, and diabetes. The secular trend of increasing obesity across all age groups has led to an increased focus on the diagnosis and management of obesity and obesity-related health risks (see Figure 1).<sup>1</sup>




Figure 1. Prevalence of obesity among adults aged 20 and over, by sex and age: United States, 2009-2010

The explosion of type 2 diabetes over the past decade is perhaps the most serious consequence of the current obesity epidemic. Nearly 2 million adult Americans were diagnosed with diabetes in 2010 alone, and the CDC estimates that 79 million Americans are pre-diabetic. Diabetes is the leading cause of kidney disease, a major cause of heart disease and stroke and the seventh leading cause of death in the United States. The financial burden of the disease is staggering; in 2007 alone total expenditures amounted to \$174 billion including medical costs and indirect costs such as disability, loss of work, and premature death from diabetes.<sup>2</sup>

There is mounting evidence that the distribution and type of excess fat may be an important prognostic indicator for disease risk. Unlike subcutaneous fat whose main function is energy storage, visceral fat cells are metabolically active and impact a wide variety of clinical risk factors including fasting glucose levels, serum triglycerides, and cholesterol.<sup>3, 4</sup>

Visceral fat occurs within the envelope formed by the abdominal muscles, principally within the greater and lesser omentum that connects the abdominal organs, and in mesenteric fat. A small amount is also found retroperitoneally.<sup>5</sup> Visceral fat is more dangerous than subcutaneous fat because visceral fat cells release proteins that contribute to inflammation, atherosclerosis, dyslipidemia, and hypertension. Visceral fat is associated with metabolic risk factors and all-cause mortality in men,<sup>6</sup> and is therefore considered a pathogenic fat depot.<sup>7</sup>




Figure 2. Computed tomography image of the abdomen at the level of the fourth lumbar vertebrae with visceral cavity outlined in white (top) and visceral fat filled in (bottom).<sup>8</sup>

## Methodology

Hologic scientists have recently developed and patented<sup>9</sup> a method for measuring visceral fat from a DXA whole body scan that is highly correlated and linearly related to visceral fat measurements by computed tomography.<sup>10</sup> This development allows clinicians and researchers to classify subjects with excess visceral fat, thereby identifying the population with the greatest obesity-related health risks. Diet, exercise, surgical, and pharmaceutical interventions can now be targeted toward high-risk individuals to maximize health benefits.

The main breakthrough in using DXA for the measurement of visceral fat was the recognition that the subcutaneous fat ring, inner abdominal muscle wall, and visceral cavity are all recognizable anatomical features on Hologic Discovery high resolution DXA whole body images. Automated software algorithms identify these anatomical features at the level of the fourth lumbar vertebrae. Using appropriate modeling, the amount of subcutaneous fat in the

abdominal region can be estimated from the subcutaneous fat on each side of the abdominal cavity. This estimate of subcutaneous fat is subtracted from the total abdominal fat in the abdominal region to yield visceral fat.

DXA visceral fat area measurements obtained with this method were highly correlated ( $r=0.93$ ) and linearly related to the visceral fat area measured by computed tomography in a study of 272 black and white women.<sup>10</sup> The low prediction error ( $SEE=16 \text{ g/cm}^2$ ) and strong functional relationship reported in this study indicates that DXA and computed tomography measurements of visceral fat are substantially equivalent. The introduction of a DXA visceral fat application confers several practical and technical advantages because DXA exams are widely available, least costly, and deliver a small fraction of the radiation dose compared with computed tomography. Furthermore, existing whole body exams, e.g. exams utilized in clinical medicine and in research studies, can be reanalyzed with the new visceral fat application as it is 100% backward compatible with all modern Hologic fan-beam technology.




Figure 3. Stylized representation of DXA projected image showing anatomical features including subcutaneous fat, abdominal wall, and visceral cavity.

### Clinical Management and Interpretation

It is important to recognize that even subjects who are normal weight and have a Body Mass Index (BMI) < 25 can have a significant accumulation of visceral fat, increasing their risk for cardiovascular disease, diabetes, and other obesity-related health risks. Likewise, some overweight or obese patients may have relatively low levels of visceral fat, normal metabolic profiles, and few or no additional risk factors. A DXA visceral fat measurement may distinguish apparently normal weight subjects with high visceral fat levels and high disease risk from metabolically normal subjects with BMIs in the overweight or obese category.

Virtually all visceral fat studies in the literature report visceral fat in units of area (cm<sup>2</sup>); only Hologic Discovery DXA scanners report visceral fat in its native area format. Cross-sectional studies have suggested that visceral fat levels exceeding 100-110 cm<sup>2</sup> adversely affect the metabolic profile in women.<sup>11</sup> Pickhardt et al. found that a visceral fat area threshold of 70 cm<sup>2</sup> yielded a sensitivity, specificity, and accuracy of 83.7%, 80.0%, and 81.3%, respectively, for metabolic syndrome in women. Visceral fat area performance was somewhat lower in men using a 125 cm<sup>2</sup> threshold.<sup>12</sup>

Nicklas found that a visceral fat area of 106 cm<sup>2</sup> was associated with an elevated risk of metabolic coronary heart disease (CHD) risk factors and a visceral fat area of 163 cm<sup>2</sup> was associated with a markedly increased risk in women. The study concluded that “these values may prove useful for defining visceral obesity and for identifying women most likely to benefit from preventive interventions.”<sup>13</sup>




Figure 4. Correlation between visceral fat area measured by computed tomography (CT) and DXA.

### Visceral Fat Diagnostic Thresholds

Diagnostic thresholds for visceral fat are becoming more firmly established as further clinical and research experience is gained. Data from such population-based studies such as NHANES will supplement the knowledge base with reference values that will provide the analytical framework for interpreting DXA visceral fat measurements. In the interim, the available literature supports visceral fat thresholds of 100 cm<sup>2</sup> for increased risk and 160 cm<sup>2</sup> for high risk.




Figure 5. Visceral fat thresholds associated with metabolic risk factors for coronary heart disease.<sup>12, 13</sup>

### Summary and Conclusion

Recent government health statistics reveal that obesity has exploded to epidemic levels. The management of clinical obesity and its related health risks is a vexing clinical problem. The reliance on antiquated measurements such as BMI, scale weight, and waist circumference further compounds the problem, as anthropomorphic measures often do not adequately assess obesity-related risks. The end result is that commonly used clinical measures of obesity misclassify individuals in terms of visceral adipose tissue and disease risk.<sup>14</sup>

Improvements to traditional measures of obesity are needed, and these refinements must be widely available and cost effective. It has been previously argued<sup>15</sup> that DXA Fat Mass Index (FMI; Fat Mass/Ht<sup>2</sup>) measurements are superior to BMI for obesity classification because FMI is a gender-specific measurement of excess fat, whereas BMI is a measure of excess weight. Likewise, DXA visceral fat measurements may allow more accurate assessment of health risks compared with non-specific measures because visceral fat is a pathogenic fat depot.

---

## References

1. Ogden CL, Carroll MD, Kit BK, Flegal KM. Prevalence of obesity in the United States, 2009–2010. NCHS data brief, no 82. Hyattsville, MD: National Center for Health Statistics. 2012. <http://www.cdc.gov/nchs/data/databriefs/db82.pdf>
2. Centers for Disease Control and Prevention. National Diabetes Fact Sheet: national estimates and general information on diabetes and prediabetes in the United States, 2011. Atlanta, GA: U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, 2011. <http://diabetes.niddk.nih.gov/dm/pubs/statistics/#fast>
3. Bergman, RN et al. Why Visceral Fat Is Bad: Mechanisms of the Metabolic Syndrome. *Obesity* (2006) 14, 16S–19S; doi: 10.1038/oby.2006.277
4. Weltman, A, et al. Impact of abdominal visceral fat, growth hormone, fitness, and insulin on lipids and lipoproteins in older adults. *Metabolism*. 2003 Jan;52(1):73-80.
5. Freedland, ES. Role of a critical visceral adipose tissue threshold (CVATT) in metabolic syndrome: implications for controlling dietary carbohydrates: a review. *Nutrition & Metabolism* 2004, 1:12
6. Kuk JL, et al. Visceral fat is an independent predictor of all-cause mortality in men. *Obesity*. 2006;14(2):336-41.
7. Fox CS, et al. Abdominal visceral and subcutaneous adipose tissue compartments: association with metabolic risk factors in the Framingham Heart Study. *Circulation*. 2007 Jul 3;116(1):39-48.
8. Yoshizumi T, Nakamura T, Yamane M, et al. Abdominal fat: standardized technique for measurement at CT. *Radiology* 1999;211:283-286.
9. T.L. Kelly, KE Wilson. Estimating visceral fat by dual-energy absorptiometry. Patent # 7,725,153 (U.S.: Hologic, Inc., 2011).
10. Micklesfield LK, Goedecke JH, Punyanitya M, Wilson KE, Kelly TL. Dual-energy x-ray performs as well as clinical computed tomography for the measurement of visceral fat. *Obesity*. 2012 May;20(5):1109-14.
11. Brochu M, Tchernof A, Turner AN, Ades PA, Poehlman ET. Is there a threshold of visceral fat loss that improves the metabolic profile in obese postmenopausal women? *Metabolism*. 2003 May;52(5):599-604.
12. Pickhardt PJ, Jee Y, O'Connor SD, Del Rio AM. Visceral Adiposity and Hepatic Steatosis at Abdominal CT: Association With the Metabolic Syndrome. *AJR Am J Roentgenol*. 2012 May;198(5):1100-7.
13. Nicklas, BJ et al. Visceral Adipose Tissue Cutoffs Associated With Metabolic Risk Factors for Coronary Heart Disease in Women. *Diabetes Care* 26:1413–1420, 2003
14. Pou, KM et al. Patterns of Abdominal Fat Distribution. *Diabetes Care*. 2009 March; 32(3): 481–485.
15. T.L. Kelly, K.E. Wilson and S.B. Heymsfield, “Dual energy X-Ray absorptiometry body composition reference values from NHANES,” *PLoS One*, 4 (2009), e7038.

### United States / Latin America

35 Crosby Drive  
Bedford, MA 01730-1401 USA  
Tel: +1.781.999.7300  
Sales: +1.781.999.7453  
Fax: +1.781.280.0668  
[www.hologic.com](http://www.hologic.com)

WP-00062 (06/12) US/International ©Hologic, Inc. 2012 All rights reserved. Printed in the USA. Specifications subject to change without notice. Hologic, Discovery and associated logos or trademarks are registered trademarks of Hologic, Inc. and/or subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshow, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to [womenshealth@hologic.com](mailto:womenshealth@hologic.com).

---